PHILIPS

RESPIRONICS

DreamWear

Gel Pillows Mask

Instructions for use - **EN**

1 Safety Information

Intended Use

This mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single patient use in the home or multi-patient use in the hospital/institutional environment. The mask is to be used on patients (>66lbs/30kg) for whom CPAP or bi-level therapy has been prescribed.

⚠ Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.


⚠ Note: This mask is not made with natural rubber latex or DEHP.

⚠ Warnings

- This mask is not intended for patients who are dependent on mechanical ventilation for their life support.
- Use of this mask while the system is not on and operating may cause rebreathing of exhaled air. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation.
- Discontinue use and contact your healthcare professional if any of the following symptoms occur: skin redness, irritation, discomfort, blurred vision, or drying of the eyes.
- Discontinue use and consult your healthcare professional if any of the following symptoms occur: unusual chest discomfort, shortness of breath, severe headache, eye pain or eye infections.
- Discontinue use and consult your healthcare professional if any of the following symptoms occur: tooth, gum, or jaw soreness. Use of a mask may aggravate an existing dental condition.

- Do not overtighten the headgear straps. Watch for signs of overtightening, such as excessive redness, sores or bulging skin around the edges of the mask. Loosen the headgear straps to alleviate symptoms.
- Do not block or seal off the exhalation ports.

2 What's in the Package

- Ⓐ Elbow with built-in exhalation (do not block)
- Ⓑ Mask frame (available in small (SM), medium (MED), large (LG))
- Ⓒ Mask pillows cushion tip
- Ⓓ Mask pillows cushion blue gel base (Do not insert into nostrils)
- Ⓔ Mask pillows cushion platform (available in S, M, L) with built-in exhalation (do not block)
- Ⓕ Mask frame fabric sleeves (optional)
- Ⓖ Headgear
- Ⓗ Tubing quick release (may remain in the CPAP hose)

3 Before Use

- Hand wash the entire mask before use (see **Mask Care** section).
- Wash your face. Do not use moisturizer or lotion on your hands or face.

- Inspect the entire mask. Discard and replace any damaged or worn parts.


4 Achieving the Right Fit

Cushion Sizing

Three cushion sizes (S, M, L) are available.

1. Place the pillows cushion tips into your nostrils.
2. Choose the cushion so the pillows cushion blue gel bases sit comfortably under the nose without any gaps. ①

Note: Do not push the pillows cushion blue gel bases into your nostrils.

(A) Mask pillows cushion blue gel base


Philips Respiration DreamStation Mask Type and System One Resistance Control Settings

Using your mask with a Philips Respiration DreamStation or System One device provides optimal comfort. The provider sets this value on your device.

Cushion Size	Setting
S, M	X2
L	X1


Putting on the Mask

1. With the mask assembled (see *Assembly* section), place the cushion under the nose and insert the pillows cushion tips into your nostrils. ②
2. Position the frame on the top of the head. ③
3. Pull the headgear over the back of the head. ④

Adjusting the Mask



1. Peel the headgear tabs away from the fabric. Adjust the straps length. Press the tabs back against the fabric to reattach. ⑤
2. Note: Do not overtighten the headgear. Signs of overtightening include redness, sores, or bulging skin around the edges of the mask.
2. Position the mask until it fits comfortably. ⑥
3. When finished the elbow should rest at the top of the head. ⑦


Mask Frame Sizing

- Medium (MED) frame
- Small (SM) frame
- Large (LG) frame

The medium (MED) mask frame will comfortably fit most faces. If the MED frame does not fit your face, contact your healthcare provider to see if a small (SM) or large (LG) mask frame would better suit your needs.

Tip: If the mask frame falls backward and is too close to your ears you may need a smaller mask frame.

Tip: If the mask frame falls forward on the head and is too close to your eyes you may need a larger mask frame.


Using the Mask

1. Connect the CPAP tubing (included with the CPAP device) to the tubing quick release on the elbow.
2. Turn the therapy device on. Lie down. Breathe normally.

3. Assume different sleeping positions. Move around until comfortable. If there are any excessive air leaks, make final adjustments. Some air leaking is normal.

Tubing Quick Release

The elbow is equipped with a tubing quick release. Pull on the tubing quick release and it will disengage from the elbow and come loose with the tubing. (8)

Removing the Mask

To keep your adjustments, remove the mask by grasping the cushion and pulling forward away from the nose. Then pull the cushion and mask up and off of the head. (9)

5 Mask Care

Cleaning Instructions

Daily: Hand wash the non-fabric parts.

Weekly: Hand wash the fabric parts.

1. Disassemble the mask (see *Disassembly* section).
2. Submerge and hand wash the mask in warm water with liquid dishwashing detergent.
3. Rinse thoroughly.
4. Lay the headgear and fabric sleeves flat or line dry. Make sure the entire mask is dry before use.

Note: Do not place the fabric parts into a clothes dryer.

Caution: Any deviation from these instructions may impact the performance of the product.

Caution: Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

Dishwashing Instructions

In addition to hand washing, the mask may be cleaned in the dishwasher once a week.

Caution: Use a mild liquid dish washing detergent only to wash the mask.

1. Remove the fabric parts. Do not wash the fabric parts in the dishwasher.
2. Wash in the top shelf of the dishwasher.

Note: Do not use the drying cycle on the dishwasher.

3. Air dry. Make sure the mask is dry before use.

Caution: Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

6 Assembly and Disassembly

Assembly

Tip: Align the triangles on the cushion and frame for proper assembly.

1. **Cushion:** Press into the mask frame until the cushion clicks into place. **10** **11**

Note: The end of the cushion and the opening in the mask frame are 'D' shaped. The cushion end should match the frame.

2. **Elbow:** Insert into the top of the mask frame.



3. **Tubing quick release:** Push the tubing quick release onto the elbow until it clicks into place.

4. **Fabric sleeves:** Wrap the sleeves around the mask frame and press the strips together.

Note: The sleeve seams should be on the outside of the frame.

5. **Headgear:** Slide the tabs into the slots on the mask frame and fold backwards.

Note: The Philips Respironics logo will be on the outside and facing up when correctly assembled.

Disassembly


1. **Headgear:** Undo the headgear tabs and pull through the slots in the mask frame.

2. **Fabric sleeves:** Undo the strips and remove from the mask frame.

3. **Tubing quick release:** Pull the tubing quick release from the elbow.

4. **Elbow:** Pull the elbow from the top of the frame.

5. **Cushion:** Pull the cushion from the mask frame.

7 Health Care Provider & Clinician Information

Multi-Patient Use


Refer to the Disinfection Guide for Professional Users to reprocess between patients in a clinical setting. Access the latest version of the Disinfection Guide at www.healthcare.philips.com or by contacting customer service at 1-800-345-6443 (USA or Canada only) or +1-724-387-4000.

Specifications

The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device.

⚠ Warning: Use outside of these specifications may result in ineffective therapy. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.

Pressure Flow Curve

Resistance

Drop in Pressure (cm H ₂ O)		at 50 SLPM	at 100 SLPM
Cushion size	S	1.3	4.8
	M	1.2	4.3
	L	1.1	3.8

Deadspace

Frame size:	SM	MED	LG
Cushion size:	S: 77 mL	S: 77 mL	S: 77 mL
	M: 80 mL	M: 80 mL	M: 80 mL
	L: 82 mL	L: 83 mL	L: 83 mL

Sound Levels

A-weighted Sound Power Level 28 dBA

A-weighted Sound Pressure Level at 1 m 20 dBA

Disposal

Dispose of in accordance with local regulations.


Storage Conditions

Temperature: -4° F to 140° F (-20° C to 60° C)

Relative Humidity: 15% to 95%, non-condensing

8 Symbols Glossary

Symbol	Title	Standard	Definition	ISO 7000 Reg. No
X1	Philips Respiration System One Resistance Control Value	Respironics symbol	Indicates the therapy device comfort setting.	N/A
X2				

Symbol	Title	Standard	Definition	ISO 7000 Reg. No
	Prescription Only	Respironics symbol	Indicates that the medical device is available by prescription only.	N/A
	Not Made with Natural Rubber Latex	Respironics symbol	Indicates that this mask is not made with natural rubber latex.	N/A
	Manufacturer	ISO 15223-1 (Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied).	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	3082
	Humidity Limitation	ISO 15223-1 (Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied).	Indicates the range of humidity to which the medical device can be safely exposed.	2620
	Temperature Limit	ISO 15223-1 (Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied).	Indicates the temperature limits to which the medical device can be safely exposed.	0632

Symbol	Title	Standard	Definition	ISO 7000 Reg. No
	Caution	ISO 15223-1 (Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied).	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	0434A
	Batch Code	ISO 15223-1 (Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied).	Indicates the manufacturer's batch code so that the batch or lot can be identified.	2492
	Catalogue Number	ISO 15223-1 (Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied).	Indicates the manufacturer's catalogue number so that the medical device can be identified.	2493

9 Limited Warranty

Respirronics, Inc. warrants that its mask systems (including mask frame and cushion) (the "Product") shall be free from defects of workmanship and materials for a period of ninety (90) days from the date of purchase (the "Warranty Period"). If the Product fails under normal conditions of use during the Warranty Period and the Product is returned to Respirronics within the Warranty Period, Respirronics will replace the Product. This warranty is nontransferable and only applies to the original owner of the Product. The foregoing replacement remedy will be the sole remedy for breach of the foregoing warranty.

This warranty does not cover damage caused by accident, misuse, abuse, negligence, alteration, failure to use or maintain the Product under conditions of normal use and in accordance with the terms of the product literature, and other defects not related to materials or workmanship.

This warranty does not apply to any Product that may have been repaired or altered by anyone other than Respirronics. Respirronics disclaims all liability for economic loss, loss of profits, overhead, or indirect, consequential, special or incidental damages which may be claimed to arise from any sale or use of the Product. Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

THIS WARRANTY IS GIVEN IN LIEU OF ALL OTHER EXPRESS WARRANTIES. IN ADDITION, ANY IMPLIED WARRANTY, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, IS SPECIFICALLY DISCLAIMED. SOME JURISDICTIONS DO NOT ALLOW DISCLAIMERS OF IMPLIED WARRANTIES, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS UNDER LAWS APPLICABLE IN YOUR SPECIFIC JURISDICTION.

To exercise your rights under this limited warranty, contact your local authorized Respirronics, Inc. dealer or Respirronics, Inc. at 1001 Murry Ridge Lane, Murrysville, Pennsylvania 15668, USA, or Respirronics Deutschland GmbH & Co KG, Gewerbestraße 17, 82211 Herrsching, Germany

Limited Warranty for Australia and New Zealand only

Note: For Australian and New Zealand customers this warranty replaces the warranty contained above.

Respirronics, Inc., a Philips Healthcare company warrants that the Products shall be free from defects of workmanship and materials and will perform in accordance with the Product specifications. 2. This warranty is valid for a period of 30 days from the date of purchase from an authorised Respirronics dealer. 3 If the Product is found to contain a defect of workmanship or materials or fails to perform in accordance with the Product specifications during the applicable warranty period, Respirronics will repair or replace, at its option, the defective material or part. 4. The customer is responsible for returning the product to an authorised Philips Respirronics dealer, and collecting the product from the authorised Philips Respirronics dealer after repair or replacement, at its own cost. Philips Respirronics is responsible only for the freight cost of transporting the product between the authorised Philips Respirronics dealer and Respirronics. Respirronics reserves the right to charge an evaluation and postage fee for any returned

Product as to which no problem is found following investigation. 5. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to materials or workmanship. 6. The warranty provided by Respiromics herein is not transferrable by the Buyer in the event of any sale or transfer of Products purchased by the Buyer from an authorised Respiromics dealer. 7. To the extent permitted by law, where the Buyer has the benefit of an implied guarantee under the Australian Consumer Law, but the Product is not of a kind ordinarily acquired for personal, domestic or household use or consumption Respiromics' liability shall be limited, at the option of Respiromics, to the replacement or repair of the Product or the supply of an equivalent Product. 8. To exercise your rights under this warranty, contact your local authorised Philips Respiromics dealer. A list of all authorised dealers is available at the following link: http://www.healthcare.philips.com/au_en/homehealth/distributors_index.wpd.

Alternatively, you can make a claim under this warranty by contacting Respiromics directly at: Philips Electronics Australia Limited, 65 Epping Road, North Ryde NSW 2113, Australia. Tel: 1300 766 488, Email: prcontact@philips.com. 9. The following statement is provided to a Buyer who is a "consumer" under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the good repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. 10. The following statement is provided to a Buyer who is a "consumer" under the Consumer Guarantees Act 1993, New Zealand: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

AUSTRALIAN SPONSOR DETAILS:

Philips Electronics Australia Ltd.
65 Epping Road, North Ryde, NSW 2113
Australia.

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© 2016 Koninklijke Philips N.V. All rights reserved.
Respironics Inc.