


Note: For proper mask sizing print out this page at 100%.
Do NOT select the "Shrink oversized pages to paper size" option.
The box in the upper right corner should measure 1 inch square.

866-414-9700
The CPAP Shop
www.TheCPAPShop.com



A large, rounded rectangular frame containing four teardrop-shaped nasal mask templates labeled P, S, M, and L from left to right. Each template has a horizontal double-headed arrow across its base with the text "cut out along lines" below it. In the bottom left corner, there are two small diagrams: one showing a person's face with a mask and dimensions W (width) and L (length), and another showing a side view of the mask with dimensions W and L. In the bottom center, the text "ComfortGel Blue nasal mask sizing gauge" is displayed above the URL "www.philips.com/comfortgelblue". In the bottom right corner, the Philips Resironics logo is shown.

ComfortGel Blue nasal mask sizing gauge
www.philips.com/comfortgelblue

PHILIPS
RESPIRONICS