

Softer, smaller, smarter...

Our second generation Full Face Mask, the *FlexiFit*TM432, takes patient comfort to another level. It incorporates a new *Soft Foam Cushion* that combines with the contoured Silicone Seal to provide an even better seal and more comfortable fit.

Experience a whole new level of comfort with Fisher & Paykel Healthcare today!

FlexiFit 432




NEW :: Soft Foam Cushion

The auto-contouring properties of the **FlexiFit™432 Full Face Mask** are enhanced by the addition of a new Soft Foam Cushion.





Full Face Delivery - *FlexiFit*[™]432



The FlexFit™ 432 Full Face Mask offers an improved seal and more comfortable fit with the inclusion of a *Soft Foam Cushion* that interacts with the Silicone Seal. This alternative to nasal or oral delivery comes fully assembled and ready for use.

UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE



FlexiFit[™] Technology Auto-contouring fit

With **Foam Cushion** FlexiFit[™] Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments. The Foam Cushion improves the sealing technology and comfort.



Glider™ Strap Anti-leak system

The sliding Glider[™] Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.



Under-Chin Design Extra stability

An Under-Chin fit provides stability by providing a constant point of reference and maintains the seal even with movement of the jaw during the night.



Efficient Bias Diffuser Quiet usage

Maintenance Free This advanced air diffuser system prevents air from blowing onto the sleeping partner and decreases noise that could potentially disrupt sleep.



Premium Frosted Silicone Seal Optimizing seal technology

This seal technology has a soft, smooth feel to slide onto the face, and its frosted texture optimizes the mask seal.



Extra stability

Complements the Under-Chin Design to maximize stability and eliminate the need for a chin strap.



Easy-Release Strap One-step removal

Simply pull the Easy-Release Cord for one-step removal of headgear. No ongoing headgear adjustments required.

CONTENTS

- Complete mask with headgear assembled and ready for use
- Three sizes, Small, Medium and Large, sold separately
- User friendly Fitting Guide

REORDER CODE: HC432

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER AND PAYKEL HEALTHCARE REPRESENTATIVE:

International PO Box 14 348, Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

USA Tel:
1800
463
3908
+1
949
453
4000
Fax:
+1
949
453
4001
Australia
Tel:
+61
3
9879
5022
Fax:
+61
3
9879
5232
Australia
Tel:
640
923
1
23
Fax:
640
924
Fax:
+1
949
453
4001
Australia
Tel:
+61
3
9879
5232
Australia
Tel:
6300
29
31
23
Fax:
644
6924
Fax:
+31
78
644
0914
China
Tel:
+86
20
3205
3486
Fax:
+86
20
3205
2132
France
Tel:
+33
1
6446
5201
Fax:
+33
1
6446</t